



Gaumard[®]
Simulators for Health Care Education



Trauma HAL[®] S3040.50

Point of Injury, Transport,
and Emergency

- Wireless and tetherless full-body trauma patient simulator
- Programmable difficult airway and surgical airway site
- Programmable spontaneous chest rise
- Trauma limbs with programmable bleeding and pressure sensors
- Rugged joints for dragging and carrying
- Real-time CPR monitoring and feedback

Features

General

- Full-body adult, 70 inches (178cm), 46 lbs. (20.9 kg)
- Rugged shoulder and hip joints allow dragging and carrying
- Splash-proof for CBRNE water spray decontamination
- Self-contained design; all operating components and blood reservoirs stored inside the body
- Wireless and tetherless control¹
- Fully operational on battery power (AC power supported)
- Interchangeable trauma/healthy lower left arm and leg
- Designed for fast scenario reset; automatic self-refilling 1.5-liter reservoir system
- Rugged for outdoor use

Neurologic

- Library of prerecorded speech responses
- Streaming voice³ (option S3040.50.300)
- Normal, miosis (constricted), and mydriasis (blown) pupil states
- Independent left/right pupil states simulate consensual and nonconsensual response

Airway

- Programmable normal and abnormal throat sounds synchronized with respiratory rate
- Realistic airway with teeth, tongue, epiglottis, and vocal cords
- Programmable difficult airway: laryngospasms
- Articulated neck and jaw for head tilt/chin lift
- Naso/orotracheal intubation: ETT, Supraglottic Airway, King LT™
- Surgical airway for cricothyrotomy, tracheotomy, and tracheal hook exercises
- Supports suctioning liquids from the airway and lungs
- Visible gastric distention with excess BVM

Breathing

- Bilateral chest rise with BVM ventilations
- Ventilation detection
- Needle decompression with audible hiss
- Programmable normal and abnormal breath sounds synchronized with respirations and ventilations
- Spontaneous chest rise and fall with normal and abnormal respiratory patterns
- Programmable unilateral chest rise
- Unilateral chest rise with right mainstem intubation
- Upper/lower anterior lung sounds

Cardiac

- Selectable normal and abnormal heart sounds synchronized with ECG and pulses
- Preprogrammed ECG rhythm library with customizable variations

- 4-lead ECG monitoring using real equipment (12-lead available option S3040.50.120)³
- Defibrillate and pace using real energy
- Effective compressions generate palpable pulses
- Real-time CPR feedback and monitoring
 - » Time to CPR
 - » Compression depth/rate
 - » Chest recoil
 - » Compression interruptions
 - » Ventilation rate
 - » Excessive ventilation
 - » Time to defibrillation
- CPR assisting tones
- Real end-tidal CO₂ capability³ (option S3040.50.078)
- Durable chest recoil design

Circulation

- Measure blood pressure with real sphygmomanometer
- Realistic Korotkoff sounds
- Bilateral blood pressure-dependent pulses
- Carotid, radial, femoral, and pedal
- Femoral artery pressure sensor
- Bilateral intravenous training arms
- Bilateral deltoid and quadriceps IM injection sites
- Intraosseous access at right tibia

Trauma

- Trauma limbs: Interchangeable trauma lower left arm and leg included
- Bleeding groin wound detects applied pressure
- Trauma site arterial bleeding synced with heart rate and blood pressure
- Blood pressure-dependent bleeding rate
- Arm and leg tourniquet detection
- Pressure on femoral artery reduces or stops distal bleeding in leg
- Internal 1.5-liter blood reservoir
- Automatic, self-refilling reservoir system
- Responsive blood loss monitor and reservoir level

Other

- Supports esophageal intubation and gastric suctioning
- Visible gastric distension with esophageal intubation

UNI® 3 Unified simulator control software

- Unified interface design
- Integrated patient monitor view
- Scenario designer
- Preprogrammed scenario library
- Real-time CPR monitoring and feedback
- Lab report designer
- Time-stamped event logging

Trauma HAL® S3040.50

S3040.50.PK 

Trauma HAL S3040.50, wireless tablet PC, "bump" case, UNI® 3 software license, battery charger, trauma/healthy lower left arms and legs, groin wound insert, soft carrying case, hands-free blood refill bag, user guide, One-Year Limited Warranty, extended warranty plans available. Patented; other patents pending.

Real CO₂ Exhalation

S3040.50.078

Adds real EtCO₂ capability. Features 10 programmable levels of CO₂ output.³

Streaming voice

S3040.50.300

Live voice streaming.³

12-Lead ECG

S3040.50.120

Adds support for 12-lead monitoring using real equipment. Includes UNI 3 Myocardial Infraction and wave editor software upgrades.³

Care in Motion™ Mobile Video-Assisted Debriefing System

CIM.PK

Care in Motion Tablet PC, 3 battery-powered HD wireless cameras, 3 adjustable camera grips, transport case, One-Year Limited Warranty, extended service plans available.

Gaumard Vitals™ Portable Virtual Monitor

30081003A

Tablet PC-based virtual patient monitor. Features 20+ dynamic scalars and waveforms, customizable alarms and layouts, and virtual defibrillator.

Request a quote

www.gaumard.com/quote
sales@gaumard.com
 Toll-Free USA & Canada
 1.800.882.6655
 Worldwide 305.971.3790